



Trafalgar Cabinet Organic Peroxide Storage 100L TCOPS100L

Self closing mechanism. Easily replaced should door be damaged

Welded shelf brackets that fully interlock with shelves.

Fully adjustable galvanised shelves perforated for any leakages and to permit free movement of air. Easily replaced if damaged

Self closing non latching door held shut by a magnetic latch which releases in the event of an internal buildup of pressure

Ventilation port with built in flame arrestor

Full instructions on usage of the cabinet attached to the inside of the door

Round corners on door.

Ventilation port with built in flame arrestor

Danger alert sign indicating specific hazardous substances stored in cabinet

Easy grip 'D' handle for opening of cabinet.

Continuous piano hinge ensuring smooth closure of door

Built in ground connector

Fully welded liquid tight sump capable of containing the entire contents of the cabinet

Painted inside and out eliminating corrosion issues.



Part Number	Capacity (L)	Shelves	External Height (mm)	External Width (mm)	External Depth (mm)	Internal Height (mm)	Internal Width (mm)	Internal Depth (mm)	Weight (kg)
BCOPS100L	100	3	1830	1110	500	1425	1015	420	191

Australian Standard

All cabinets comply with AS2714- 1993 as follows.

1. Each cabinet is marked with the name and address of the manufacturer
2. The maximum storage capacity.
3. A sign bearing the words "NO SMOKING, NO IGNITION SOURCES WITHIN 3M" in lettering at least 50mm high

Storage

Suitable for safe storage of organic peroxide materials as defined in the ADG code.

Options

Full range of spare parts including additional shelves, closing mechanisms etc. available.



TRAFALGAR



CABINETS

Phone: 1800 888 714
 Email: sales@tfire.com.au
 Website: www.tfire.com.au

Trafalgar reserves the right to change specifications without notice. Please check with your supplier at the time of order.